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## Amino Graft Information and Instructions for Use

### Description

Amnio Graft is a Syrian product prepared and controlled by quality system established and documented in quality assurance office in Syrian Atomic Energy Commission (SAEC).

Amnio Graft is a wound, burn and ulcer dressing including dried human sterilized amniotic membrane, contains a human tissues, cells, and cellular product prepared from placentas obtained from Cesarean section, and prepared according to current Good Tissue Practices (cGTP), Good Manufacturing Practices (cGMP), and Good Processing Practices (cGPP).

### Amnio Graft preparation

- The amniotic membrane was obtained from healthy donor determined to be eligible based on the results of screening and testing to determine if any risk factor has been found, and documented by maternity hospital and Syrian Atomic Energy Commission laboratories. All infectious disease tests should be negative to be the amniotic membrane suitable for preparation the Amnio Graft.
- The Amnio Graft was prepared in amniotic membrane facility production in radiation technology department, Syrian Atomic Energy Commission (SAEC) according to the regulations recommended by National Committee and established by the Syrian Atomic Energy Commission (SAEC).
- The final product of Amnio Graft was sterilized by gamma irradiation at dose ranged from 25 to 35 kGy according to the code of practice for the radiation sterilization of tissue allograft. That recommended as sterilization dose from International Atomic Energy Agency (IAEA) and World Health Organization (WHO).

### Indications for use (Clinical Applications)

Amnio Graft is an allograft intended for use in the management of interior or exterior acute and chronic, wounds as well as appropriate surgical applications. including

wounds, Burns, Diabetic foot ulcers (DFUs), periodontal disease, cancer, and osteoarthritis. All infectious disease tests were negative. This allograft tissue has been determined to be suitable for transplantation.

### Amnio Graft components properties

- Contains different types Collagen, Laminin, Elastin Fibronectin, Kproteoglycans, Nidogens, Perlecan, Agrin, and Cytokines.
- Contains Non-immunogenic components.
- Stem cells source.
- Anti-inflammatory, antimicrobial and antiviral properties.
- Secretion of growth factors and lack of anti-genicity.
- Amnio Graft is flexibility available Applicable Detachable without Trauma, Biocompatible, Economic, Nontoxic, and Non-allergenic.

### Amnio Graft benefits

- Acceleration of wound, burn and ulcer healing.
- Cellular occlusive with rapid sealing of protected space.
- Shown to recruit mesenchymal stem cells
- Aids in the formation of clinical attachment.
- Allows for rapid establishment of blood supply.
- Protect the injury area from external infection and suppresses inflammation.
- Anti-scarring and anti inflammatory properties.
- Minimizing fluid/water loss (Prevent dehydration).
- Reduces losses proteins.
- Extra exudates.
- Pain reduction (Alleviation of pain).

### Storage

Amnio Graft is a dried and sterilized amniotic membrane could be stored at ambient temperature. Do not refrigerate or freeze. It is the responsibility of the transplant facility or clinician to maintain the tissue intended for transplantation in the appropriate recommended storage conditions prior to transplant. Within the expiration date printed on

product packaging (shelf-life is 2 years from date of manufacture).

### Instructions

- Open the pouch to retrieve the sterilized packaged Amnio graft.
- Open the clear inner peel pouch and retrieve Amnio Graf using aseptic techniques.
- Aseptically transfer approximately 0.5-2.0 mL of normal saline or sterile injection water into the amniotic membrane with a sterilized needle.
- Remove the tissue from the inner pouch using sterile gloves/forceps.
- Place the tissue on the surgical area to deliver the therapeutic actions of AmnioGraft.

**Note:** Amnio Graft is for single use only in one patient by a licensed physician (DPM, MD, or DO) or a qualified medical professional with relevant clinical experience such as a nurse or PA.

### Precautions

- Do not use Amnio Graft if the internal package has been violated, opened, or damaged, or if mishandling has caused possible damage or contamination.
- Do not use Amnio Graft if there is any abnormality observed (e.g. labeling, shipping, missing information, etc.).
- Do not use if expiration date has been exceeded.
- Do not autoclave the product for sterilization before use.
- Do not sterilize or re-sterilize by any methods.
- The healthcare professional is responsible for informing the patient of the risks associated with his/her treatment and the possibility of complications or adverse reactions.

**Note:** No antibiotics were used during the processing of Amnio Graft.

Amnio Graft produced and sterilized in Human Amniotic Graft Producing Facility, Radiation Technology Department, Syrian Atomic Energy Commission. Damascus, P.O Box 609.

**Fax:** +963116113389  
**Phone:** +963112132580  
**E-mail:** [atomic@aec.org.sy](mailto:atomic@aec.org.sy)